




SENSUS

An at-home screening device for diabetics at risk of Peripheral Neuropathy (PN)




Early detection of PN is essential for disease mitigation

Symptoms of PN

-  Numbness
-  Imbalance
-  Sensitivity

Purchase Process

- 1 Doctor recommends Sensus
- 2 Manufacturing & customization
- 3 Guided user testing
- 4 At-home personal use



Specifications

Weight: 4.75 lbs
Size: 14" x 14" x 3.7"
Force Applied: 0.5 - 0.7 N

13%

Of adults have diabetes in the USA

50%

Of these diabetics will have some form of PN

